

Document Title: <b>Slick, CheckLite 250 Plus</b>		Part # and Rev. <b>18082</b>	
		Release Date:	
Rev.	Description	Revision By	Date
Date -- Time	Added description from other slicks	Allison	6/9/16

**Production Filename:**

[http://publications.decagon.com/Marketing/Description%20Files/18082\\_Slick\\_CheckLite\\_250\\_Plus.pdf](http://publications.decagon.com/Marketing/Description%20Files/18082_Slick_CheckLite_250_Plus.pdf)

**Printing Process:** Digital

**Size:** 8.5 inches wide, 11 inches tall

**Ink:** 4/4 Color Both Sides

**Paper:** 100 lbs. book gloss coated

**Page Count:** 2

**Proof type:** soft proof

**Single hole drilling:** none



Kit for microbial biomass assays

## CheckLite™ 250 Plus



### Characteristics of kit

- CheckLite™ 250 Plus is a kit for microbial biomass assays based on the ATP measuring method.
- CheckLite™ contains thermostable firefly luciferase.
- CheckLite™ is characterized by its ease of use and rapid performance.
- CheckLite™ can detect 1,000 cells of coliforms/ml (100 cell/assay).

### Composition of kit

(250 assays/kit)

Luciferin-luciferase reagent (lyophilized)

5

Reconstitution buffer (5.5 ml)

5

ATP releasing reagent (5.5 ml)

5

CheckLite™ 250 Plus should be stored at 2-8°C in a refrigerator. DO NOT FREEZE.

Manufactured by

 **Kikkoman Corporation**

250 Noda, Noda-city Chiba-pref., 278-8001 JAPAN

Inquiries

Biochemicals Division

Tel: +81-3-5521-5490 Fax: +81-3-5521-5498

Email: Biochem@mail.kikkoman.co.jp

Specifications are subjected to change without notice.

©2001 Kikkoman Corp. (06-1)

