

Document Title: Slick, CheckLite HS Set		Part # and Rev. 18083	
		Release Date:	
Rev.	Description	Revision By	Date
Date -- Time	Added description from slick	Allison	6/8/16

Production Filename:

http://publications.decagon.com/Marketing/Description%20Files/18083_Slick_CheckLite_HS_Set.pdf

Printing Process: Digital

Size: 8.5 inches wide, 11 inches tall

Ink: 4/4 Color Both Sides

Paper: 100 lbs. book gloss coated

Page Count: 2

Proof type: soft proof

Single hole drilling: none

kikkoman [®] Kit for highly sensitive microbial biomass assays

CheckLite™ HS Set

Luciferin-luciferase ATP assay reagent kit for hypersensitive viable count assay.



Characteristics of kit

- CheckLite™ HS Set is a kit for microbial biomass assays based on the ATP measuring method.
- CheckLite™ contains thermostable firefly luciferase and ATP eliminating enzymes.
- CheckLite™ is characterized by its ease of use, rapid performance and high sensitivity.
- CheckLite™ can detect 100 cells of coliforms/ml (10 cell/assay).

Composition of kit	(100 assays/kit)
Luciferin-luciferase reagent HS (lyophilized)	2
Reconstitution buffer (5.5 ml)	2
ATP eliminating reagent (lyophilized)	2
Reconstitution buffer (5.5 ml)	2
ATP releasing reagent (5.5 ml)	2

Example of biomass assays



CheckLite™ HS Set should be stored at 2-8°C in a refrigerator. DO NOT FREEZE.