

# CheckLite™ 250 Plus



## Characteristics of kit

- CheckLite™ 250 Plus is a kit for microbial biomass assays based on the ATP measuring method.
- CheckLite™ contains thermostable firefly luciferase.
- CheckLite™ is characterized by its ease of use and rapid performance.
- CheckLite™ can detect 1,000 cells of coliforms/ml (100 cell/assay).

## Composition of kit

(250 assays/kit)

Luciferin-luciferase reagent (lyophilized)	5
Reconstitution buffer (5.5 ml)	5
ATP releasing reagent (5.5 ml)	5

CheckLite™ 250 Plus should be stored at 2-8°C in a refrigerator. DO NOT FREEZE.

Manufactured by

**KIKKOMAN** **Kikkoman Corporation**  
250 Noda, Noda-city Chiba-pref., 278-8061 JAPAN  
Inquiries

**Biochemicals Division**  
Tel. +81-3-5521-5490 Fax +81-3-5521-5498  
Email: Biochem@mail.kikkoman.co.jp

Specifications are subjected to change without notice.  
©2001 Kikkoman Corp. (06-1)

