



METER

September 10, 2026

Re: Sole Source Certificate for WP4C Dewpoint Potentiometer

To Whom It May Concern:

This is to certify that METER Group, Inc., a company based in Pullman, WA, is the sole manufacturer of the WP4C Dewpoint Potentiometer.

The WP4C uses the chilled mirror dew point method to determine water potential automatically. This makes the WP4C an accurate and fast water potential instrument option. Other instruments measure relative humidity with sensors that detect changes in electrical resistance or capacitance; still others, like pressure plates and hanging water columns, use outdated manual pressure measurements, and are not accurate in the range of the WP4C. The WP4C is compatible with LABROS SoilView software and can be used with HYPROP data to generate a complete soil moisture retention curve.

The WP4C is manufactured in Pullman, Washington. In the United States, it is currently sold through METER Group; in other markets, it is available through METER Group's established distribution channels. Each instrument's output is verified at our facility before it is shipped. The product has a one-year warranty against manufacturing defects. It is recommended to have your WP4C unit cleaned and verified annually by METER Group.

Any comments or questions regarding this document can be brought to my attention.

Sincerely,

Shaun Weldon

Product Scientist